Introducing new implant technology



- New wide Ponto implant
- New 12 mm Ponto abutment

NEW

New ISQ compatible connection screw

New wide Ponto implant - designed to improve stability

OptiGrip implant geometry

The industry's largest initial implant surface in contact with bone •.

Universal hexagon interface

Full backwards compatibility

Proven Brånemark® type surface

Wider (Ø 4.5 mm) diameter

New 12 mm Ponto abutment

Improved adaptability to different skin thicknesses and new surgical techniques





Green ISQ compatible connection screw

Easy to identify and allows for monitoring of implant stability

